

Full custom configurations

Tomtec, Inc. is a full service US manufacturer with its in-house CAD, tool & die design manufacturing and engineering support. This allows us to quickly design, develop and deploy custom manufactured components and/or fixtures to meet your demanding work-flow needs.

	Pipette Head Configuration	Quadra4 hardware Options	Station Accessories
High Throughput Screening for NME	 96 channel DMV pipette head 384 channel 60µL Stainless steel pipette head 	 Quadra indexing stage for Serial dilution & Reformatting with Individual Tip Loader Bi Directional Stacker with temp. control for continuous runs Enhanced Pneumatic Control System N₇ Blow-through with vacuum control Thermal lamination unit to seal assays Sample Identification Module Barcode Reader for LIMS integration Automated vacuum filtration 	 Temperature controlled heating station Shaker Nest Assembly Constant Level Re-circulating Reservoir System - Includes: Reservoir, Peristaltic Pump, and Tubing Kit.
Solid Phase Extraction	 96 channel 450μL pipette head 96 channel DMV pipette head 	 Quadra indexing stage for Serial dilution & Reformatting with Individual Tip Loader Enhanced Pneumatic Control System Auto Regulated Vacuum Trap Automated Filtration System Sample Identification Module Barcode Reader for LIMS integration Bi Directional Stacker for continuous runs 	 Ultrasonic tip wash System Shaker Nest Assembly Heat/Chill Nest (Peltier) Positive Pressure SPE
Sandwich ELISA	96 channel 450μL pipette head	 Quadra indexing stage for Serial dilution & Reformatting with Individual Tip Loader Enhanced Pneumatic Control System Sample Identification Module Barcode Reader for LIMS integration 	 Ultrasonic tip wash System Shaker Nest Assembly 3 - Heat/Chill Nest (Peltier)
Genomics - sample preparation and PCR	96 channel 450μL pipette head	 Enhanced Pneumatic Control System Sample Identification Module Barcode Reader for LIMS integration Auto-Mag (Magnetic Plate Nest) 96 Well 	 3 stainless steel Coated Aluminum Reservoir Ultrasonic tip wash System Shaker Nest Assembly Heat/Chill Nest (Peltier)
Drug testing (metabolites - NIDA-5)	• 96 channel 60µL stainless steel pipette head	 Quadra indexing stage for Serial dilution & Reformatting with Individual Tip Loader Enhanced Pneumatic Control System Sample Identification Module Barcode Reader for LIMS integration 	 Ultrasonic tip wash System Constant Level Re-circulating Reservoir System - Includes: Reservoir, Peristaltic Pump, and Tubing Kit. Shaker Nest Assembly