

Pipette Replacement Tips

Clear cannula Polypropylene Pipette Tips for the Quadra & Formatter Family of Liquid Handling workstations

Tip Description	P/N	Quadra 450uL Head	Quadra DMV ¹ Head	Quadra DSV ² Head	Formatter & Formatter2
450µl Capacity, Standard Orifice (.036) 96 tips/rack 50 racks/carton Non-Sterile *	196-205	X			X
450µl Capacity, Standard Orifice (.036) 96 tips/rack, 50 racks/carton Sterilized *	196-205S	X			X
20µl Capacity 96 tips/rack, 5 racks/package, 10 packages/carton Non-Sterile	196-305			X	
20µl Capacity 96 tips/rack 5 racks/package, 10 packages/carton Sterile	196-305S			X	
250µl Capacity 96 tips/rack, 10 racks per Box, 5 Boxes per Carton Non-Sterile*	196-306		X		
50µl Capacity 96 tips/rack 10 racks per Box, 5 Boxes per Carton. Non-Sterile*	196-307		X		
20µl Capacity 96 tips/rack 10 racks per Box, 5 Boxes per Carton Non-Sterile *	196-308		X		

* Indicates volume discounts are available, please call or email for quote

(1) Disposable Multi-Volume Pipette Tips

(2) Disposable Small-Volume Pipette Tips